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Operate Manual

目录 CATALOGUE

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The screenshot shows a web interface for login and registration. At the top left, there is a 'Language' dropdown menu set to '简体中文'. At the top right, there is a 'Wifi设置' button. Below the header, there are two tabs: '登录' (Login) and '注册' (Register), with '注册' being the active tab. The registration form consists of several input fields and buttons:

- Input field: 请输入账号 (Please enter account)
- Input field: 请输入公司名称 (Please enter company name)
- Input field: 请输入密码 (Please enter password)
- Input field: 请输入公司地址 (Please enter company address)
- Input field: 请确定密码 (Please confirm password)
- Input field: 请输入联系人姓名 (Please enter contact name)
- Input field: 请输入手机号码 (Please enter mobile number)
- Input field: 请输入邮箱 (Please enter email)
- Input field: 请输入验证码 (Please enter verification code)
- Button: 获取验证码 (Get verification code)
- Button: 注册 (Register)

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The screenshot shows a software interface for calibration. At the top, there is a navigation bar with a home icon, the word '校准' (Calibration), and 'ASTM C/2°'. Below the navigation bar, there is a dropdown menu for '校准有效时间' (Calibration valid time) set to '8小时' (8 hours). The main area is divided into two sections. On the left, there is an image of a haze meter device. On the right, the text reads '第二步 100%' (Step 2 100%), followed by instructions: '请保持haze-port对准空气，并点击测量键进行校准' (Please keep the haze-port aligned with the air and click the measurement key for calibration). Below the instructions, it says '仪器0%校准成功' (Instrument 0% calibration successful) in green. At the bottom right of the main area, there is a grey button labeled '校准' (Calibration).

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The screenshot displays the '个人中心' (Personal Center) interface. At the top, there is a navigation bar with a home icon, the text '个人中心', and the user's name 'ASTM C/2°'. Below the navigation bar, the interface is divided into two main sections. The left section, titled '账号: admin', contains four input fields for '公司名称' (Company Name), '公司地址' (Company Address), '联系人姓名' (Contact Name), and '邮箱' (Email). The right section, titled '修改账户密码' (Change Account Password), contains three input fields for '原密码' (Original Password), '新密码' (New Password), and '确认密码' (Confirm Password), along with a '确定' (Confirm) button. At the bottom right of the interface, there is a '注销' (Logout) button.

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
The screenshot shows a software interface for color measurement. At the top right, it says 'ASTM C/2°'. Below that, there are two tabs: '标样测量 -- 标样0003' and '试样测量 -- 试样0001'. The interface is divided into three main sections:

- 标样 (Standard):** L* = 82.03, a* = 1.46, b* = 7.30, c* = 7.45, h = 78.66
- 试样 (Sample):** L* = 82.00, a* = 1.54, b* = 7.25, c* = 7.41, h = 78.04
- Comparison Results:**
 - dL* = -0.04 合格
 - da* = 0.07 合格
 - db* = -0.05 合格
 - dc* = -0.03 合格
 - dH* = -0.08 合格
 - dEcmc(2.0:1.0) 0.13 合格

On the right side, there is a vertical toolbar with icons for '设置' (Settings), '保存' (Save), and '测量' (Measure).

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The screenshot displays the Hunter Lab software interface. At the top, it shows 'ASTM C/2°'. Below this, there are two tabs: '标样测量 --- 标样0003' (Standard Measurement --- Standard 0003) and '试样测量 --- 试样0001' (Sample Measurement --- Sample 0001). The interface is divided into three main sections:

- 标样 (Standard):** (Hunter) L = 77.67, a = 1.48, b = 6.70
- 试样 (Sample):** (Hunter) L = 77.66, a = 1.48, b = 6.68
- 差值 (Differences):** dL = -0.01 (合格), da = -0.01 (合格), db = -0.01 (合格). dEab = 0.02 (合格).

On the right side, there is a vertical menu with icons and labels: '设置' (Settings), '保存' (Save), and '测量' (Measure).

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




The screenshot displays the R-6 software interface for absorbance measurement. It features a graph on the left and a data table on the right. The graph, titled '吸光度曲线', plots absorbance (y-axis, 0 to 5) against wavelength in nanometers (x-axis, 400 to 700). Two curves are shown: a green line for '标样' (reference) and a red line for '试样' (sample). The data table compares the absorbance values for '标样' and '试样' at various wavelengths, showing a difference of 0.00 for most points and -0.00 for 440nm and 460nm. The interface also includes a '设置' (Settings) button, a '保存' (Save) button, and a '测量' (Measure) button.

波长	标样	试样	差值
400nm	0.00	0.00	0.00
410nm	0.00	0.00	0.00
420nm	0.00	0.00	0.00
430nm	0.00	0.00	0.00
440nm	0.31	0.31	-0.00
450nm	0.30	0.30	0.00
460nm	0.29	0.29	-0.00
470nm	0.29	0.29	0.00

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


Operation Manual

CATALOGUE

[I] Usage Notice-----	01
[II] Precautions-----	02
[III] Functional description-----	03
[IV] Technical Paramet -----	04
[V] Introduction to appearance structure-----	05
[VI] Measurement flow chart-----	06
[VII] Software interface introductio -----	08
7.1 Features-----	08
7.2 Login-----	10
7.3 Guide-----	13
7.4 Calibration-----	19
7.5 Measure-----	22
7.6 Set up -----	25
7.7 Data browsing -----	32
7.8 Personal center-----	37
7.9 About -----	38
7.10 Log -----	39
7.11 Update-----	40
[VIII] Introduction to the measurement interface-----	41
8.1 Haze/transmittance/Clarity -----	42
8.2 Haze/transmittance -----	42
8.3 Haze-----	43
8.4 Transmittance -----	44
8.5 Clarity -----	45
8.6 Color difference-----	46
8.7 Data -----	52
8.8 Figure -----	53
8.9 Metamerism-----	59
8.10 Liquid-----	60
[IX] Parameter introduction-----	61
[X] Exception handling analysis-----	65
[XI] Annex -----	66
[XII] Company statement-----	68

Usage Notice

- 1、 This machine is an instrument designed to measure the haze, total light transmittance, and spectral transmittance of transparent, translucent parallel plane materials such as plastics, films, glass products, LCD panels, etc.
- 2、 This machine is widely used in laboratory, factory, or field operation, enough to achieve the best haze and color measurement in quality control in almost all application fields.
- 3、 The time period of the limited warranty is the time since the purchase of this instrument (time: such as one year). If your instrument needs service, please bring the instrument to the local sales department to contact us for repair.
In order to avoid affecting the accuracy of the instrument, please do not
- 4、 Disassemble the instrument without permission. If the instrument is damaged due to unauthorized disassembly of the machine or incorrect use, the user is responsible for it.

Precautions

- 1、 Put the instrument on the flat place carefully.
- 2、 The instrument is not moistureproof, please keep the instrument in the dry area.
- 3、 Large force, or sharp objects may damage the screen.
- 4、 It is recommended to use the original power adapter which comes with the instrument.
- 5、 To ensure the machine to work properly, please do not store, or use the instrument in places that are too hot or too cold; please do not put the machine in damp locations, or directly under sunlight. Do not use the instrument in severe environment such as strong shock or quake.
- 6、 Please avoid strong electromagnetic interference in usage.
- 7、 To ensure accurate measurement, please keep the instrument stable during the test and do not shake it.
- 8、 Instrument belongs to high precision instrument, please keep it well after usage.
- 9、 Do not clean the inside of the integrating sphere.
- 11、 If the instrument fails, please don't try to repair it by yourself, our customer service department will quickly help the customer.

12、 Any future update on the manual, we are not obliged to notify you.
If any questions, please contact us directly.



Functional description

- 1、 Color and Haze meter conforms to standards: GB/T 2410, ASTM D 1003, ISO 13468, ISO 14782 and GB/T 2410-2008, JJF 1303-2011.
- 2、 Color and Haze meter can provide haze and transmittance value under CIE-A, CIE-C and CIE-D65 light sources.
- 3、 Color and Haze meter contains compensation method which provides higher accuracy on test result.
- 4、 Open sample measurement area to make it can measure samples at any sizes.
- 5、 Haze meter adopts 7.0 inch capacitive touch screen and Android operating system for friendly interface.
- 6、 Haze, color and transmittance analysis software to meet customer's requirement on data analysis and management.

Technical Parameters



Standard	ASTM D1003/D1044,ISO13468/ISO14782,JIS K 7105,JIS K 7361,JIS K 7136,GB/T 2410-08,CIE No.15,ISO 7724/1,ASTM E1164,DIN 5033 Teil7,JIS Z8722 Condition c
Geometry	0/d
Light Source	Haze/ Transmittance: CIE-A,CIE-C,CIE-D65 Chromaticity Value: A,C,D50,D55,D65,D75,F1,F2,F3,F4,F5,F6,F7,F8,F9,F10,F11,F12,CMF, U30,DLF,NBF,TL83,TL84
Spectral Response	CIE Luminosity function Y/V(λ)
Measurement Parameter	HAZE,Transmittance (T),CIE Lab,LCh,CIE Luv,XYZ,Yxy,Spectral Transmittance,Clarity, Hunter Lab MunsellMI,CMYK,WI(ASTM E313-00,ASTM E313-73,CIE/ISO,Hunter, Taube Berger Stensby),YI(ASTM D1925,ASTM E313-00,ASTM E313-73), Tint(ASTM E313-00),Metamerism IndexMilm,APHA,Pt-Co,Gardner, Color difference (ΔE*ab, ΔE*CH, ΔE*uv, ΔE*cmc(2:1), ΔE*cmc(1:1), ΔE*94, ΔE*00)
Wavelength	400-700nm
Wavelength Interval	10nm
Measurement Range	0-100%
Haze Resolution	0.01unit
Haze Repeatability	≤0.1unit
Measurement Area/Sample Port	18mm/25.4mm
Sample Size	Thickness≤145mm
Display	7 inch capacitive touch screen
Memory	20000values
Interface	USB
Power	220V (with power cable)
Working Temperature	0~45℃, Relative humidity 80% or lower (35℃) no condensation
Storage Temperature	-25℃~55℃, Relative humidity 80% or lower (35℃) no condensation
Size	LxWxH: 598mmX247mmX366mm
Standard Accessory	PC software (Haze QC)
Optional	Measurement Fixture、Haze Standard Plate、5mm/7mm/10mmCaliber plate、40*10Cuvette

Appearance and structure


- ① Display Screen ② Light Exit Aperture ③ Compensation Aperture ④ Test Aperture
⑤ USB ⑥ Switch on-off ⑦ Power Interface ⑧ Light Source Exchange port


Measurement flow chart


0% Calibration, 100% Calibration, Clarity Calibration

Target Measure and Sample Measure

Check measurement results

Save/not save the test result

Finish


Software interface introduction

[Features]

A-1

Main Interface


The software is divided into 8 modules, namely: calibration, measurement, setting, data browsing, personal center, about, log, update.

A-2

Title

From left to right, the title bar is: home button, measurement mode, light source and viewing angle.

[Log In]

Login is divided into local login and network login. Check the Remember password, the account and password will be automatically entered when the next boot, check from The next time you start automatic login, you will skip the login page and go directly to the software.

B-1

The local login account is: admin, and the password is the instrument serial number by default (you can modify it in the personal center after login). For example, the instrument serial number is: C81118C0128, then enter C81118C0128 at the login password.

The screenshot shows a login page with a dark header. On the left, there is a 'Language' dropdown menu set to 'English'. On the right, there is a 'Wifi Setting' button. Below the header, there are two tabs: 'Login' (which is active) and 'Register'. The main content area contains a 'Local Login' dropdown menu. Below it are two input fields: the first contains the text 'admin', and the second contains a series of dots representing a masked password. There are two checkboxes: 'Keep Password' (checked) and 'Auto Login' (unchecked). A blue 'Login' button is positioned below the checkboxes. At the bottom right of the form area, there is a link that says 'Forget Password?'.

B-1

B-2

Network Login

You need to connect to the Internet to log in to the Internet. Click Register to register your account. You can choose your email or mobile phone number to register. After the registration is completed, you can log in using the registered account. Use the network login to upload the stored data to the cloud, and perform data management on the windows.


The screenshot shows a web interface with a dark header. On the left, there is a 'Language' dropdown menu set to 'English'. On the right, there is a 'Wifi Setting' button. Below the header, there are two tabs: 'Login' and 'Register'. The 'Register' tab is active. The registration form consists of several input fields: 'Account', 'Company Name', 'Password', 'Address', 'Confirm Password', 'Name', 'Cell Phone', and 'Email'. At the bottom of the form, there are three buttons: 'Input Code', 'Get Code', and 'Register'.

B-2

[Guide]

C

After logging into the software for the first time or re-login after restoring the factory, you will enter the instrument user guide, and you can set the instrument according to the guide.

Instrument Setting

Exit

Mode Selection

ASTM

ISO

System Setting

Screen Light

Language

English

Factory Reset

Screen Rotation

Wifi Setting

Timezone Setting

1/5

Previous

Next

C-2

Parameter Setting

Exit

Illuminant&Angle (The second illuminant is for calculating the metamerism)

First

C

2°

Second

A

2°

CMC(l:c)

l 2.0

c 1.0

CIE94

KL

1.0

KC

1.0

KH

1.0

CIE 2000

KL 1.0

KC 1.0

KH 1.0

2/5

Previous

Next

C-3

Tolerance Setting

Exit

Haze/T

CIE LAB&LCH

Hunter Lab

CIEDE2000

CIE LUV

CMC(l:c)&CIE94

Liquid

Temperature&Humidity

Haze/T

		Greater than	Less than	Between
dh(Haze) : ±	<input type="text" value="2.0"/>	<input type="text" value="Fail"/>	<input type="text" value="Fail"/>	<input type="text" value="Pass"/>
dt(T) : ±	<input type="text" value="2.0"/>	<input type="text" value="Fail"/>	<input type="text" value="Fail"/>	<input type="text" value="Pass"/>
dC : ±	<input type="text" value="2.0"/>	<input type="text" value="Fail"/>	<input type="text" value="Fail"/>	<input type="text" value="Pass"/>

3/5

Previous

Next

C-4

Other Settings

Exit

Average

Single Test Average Test

Save Setting

Manual Save Auto Save

Naming Rules

Target + Number + Date

Sample + Number + Date

4/5

Previous

Next

C-5

Display Setting

- H/T/C
- Haze/T**
- Haze
- Transmittance
- Clarity
- Color Diff.
- Data
- Figure
- Metamerism
- Liquid

5/5

Previous

Finish

C-6


[Calibration]

D-1

0% calibration

Perform 0% calibration according to the software prompt. In the ASTM mode, please cover the compensation port (comp-port) with the compensation port cap; in ISO mode, please keep the compensation port (comp-port) open.

This page can be used to set the calibration time. The default calibration valid time is 8 hours.

D-1

100% calibration

Perform 100% calibration as prompted by the software.

Calibrate

ASTM C/2°

Valid Time 8 Hours

Step Two

100%

Keep haze-port aligned with the air and click the measurement button for calibration.

0% calibration succeed

Calibrate

Clarity calibration

Perform clarity calibration as prompted by the software.

The screenshot shows a software interface for a calibration process. At the top, there is a dark grey header with a home icon, the word 'Calibrate', and 'ASTM C/2°'. Below the header, on the left, is a 'Valid Time' dropdown menu currently set to '8 Hours'. In the center, there is a 3D rendering of a white and black device with a blue touchscreen. To the right of the device, the text reads 'Step Three' followed by 'Clarity Adjust'. Below this, instructions state: 'Please place the calibration sheet on the clarity-port and click the measurement button to calibrate the clarity'. A green message '100% calibration succeed' is displayed. At the bottom right of the instruction area is a grey 'Calibrate' button.

[Measure]

The measurement is divided into three modes: standard sample measurement, sample measurement and other measurement modes (metamorphism, liquid color)

E-1

Standard measurement

Click Measure on the homepage to enter the standard sample measurement interface. After the sample is placed, click the measurement button in the lower right corner of the screen. The measurement data is displayed on the screen. The measurement button returns to the pressable state, indicating that the measurement is complete.

ASTM C/2°

Test Target --- Target0001

Test Sample

Setting

Save

Measure

Click to measure

H	T
0.00	100.05

E-1

E-2

Sample measurement

In the standard sample measurement interface, click the sample measurement on the interface to switch to the sample measurement interface. Similarly, after the sample is placed, click the measurement button in the lower right corner of the screen to perform sample data measurement.

The screenshot displays the sample measurement interface. At the top, it shows 'Test Target --- Target0001' and 'Test Sample --- Sample0001'. The main area features a 'Transmittance Curve(%)' graph with 'Wavelength (nm)' on the x-axis (400 to 700) and transmittance on the y-axis (0 to 200). A red line represents the 'Sample' and a green line represents the 'Target'. Below the graph is a color calibration bar. To the right of the graph is a data table:

	H	T
Target	0.00	100.05
Sample	0.00	100.07
Difference	0.00	0.02
Dudge	Pass	Pass

On the right side of the interface, there are buttons for 'Setting', 'Save', and 'Measure'. The 'Measure' button is highlighted with a red box and a red arrow pointing to it with the text 'Click to measure'. A blue oval with 'E-2' is also present in the bottom right corner of the interface.

[Setup]


The setting interface can configure the instrument measurement method, data calculation parameters, tolerance, software display, saving method, naming rules, average, etc. After modification, you need to click the "Apply" button.

F-1

Instrument settings

The instrument settings are divided into 6 areas:

1. Haze measurement mode: ASTM mode or ISO mode can be set;
2. System settings: You can set the screen backlight and language switch;
3. Factory reset: software configuration is restored to the initial state;
4. Screen rotation: Click on the screen to rotate the screen and the color will reverse 180°;
5. Wifi settings: you can choose WiFi and log in;
6. Time zone setting: time display in different countries, network can automatically synchronize time.

F-2

Parameter settings

Parameter setting can configure the data calculation

1. Light source & angle: You can set the light source and angle of the calculated data. The first light source and angle are the calculated data in all modes, and the second light source is only used to calculate the metamerism (Note: the data of the similar color search display is fixed at D65/ 10°)

- 2、CMC(l:c):L:c coefficient of CMC color difference formula can be set;
- 3、CIE 94:Can set KL, KC, KH coefficients of CIE94 color difference formula;
- 4、CIE 2000:Can set KL, KC, KH coefficients of CIE2000 color difference formula;

The screenshot shows a software settings window titled "Setting" with a home icon on the left and "ASTM C/2°" on the right. Below the title bar are five tabs: "Instrument", "Parameter", "Tolerance", "Display", and "Other". To the right of these tabs are two buttons: "Guide" and "Apply".

The "Parameter" tab is active and contains four main sections:


- Illuminant&Angle** (The second illuminant is for calculating the metamerism):
 - First: C (dropdown), 2° (dropdown)
 - Second: A (dropdown), 2° (dropdown)
- CMC(l:c)**:
 - l: 2.0 (input field)
 - c: 1.0 (input field)
- CIE94**:
 - KL: 1.0 (input field)
 - KC: 1.0 (input field)
 - KH: 1.0 (input field)
- CIE 2000**:
 - KL: 1.0 (input field)
 - KC: 1.0 (input field)
 - KH: 1.0 (input field)

F-2

F-3

Tolerance settings

Tolerance is used to judge whether the measurement data is qualified. When the measurement data exceeds the tolerance range, the data will be disqualified. When the measurement data is less than or equal to the tolerance, the data will be prompted to be qualified. In this interface, you can set the tolerance of different color difference formulas and modes. (Where CIE LAB can customize the prompt language)

F-3

Display setting

Display settings can set the content displayed under "Measurement Page". Divided into the following:

1. Haze/transmittance/sharpness: display haze, total transmittance and sharpness parameters at the same time;
2. Haze/transmittance: Simultaneously display the parameters of haze and total transmittance;
3. Haze: Only display haze parameters;
4. Transmittance: Only the total transmittance parameter is displayed;
5. Clarity: only display the clarity parameter;
5. Color difference: CIELABCH, CIEDE2000, CIE94, CMC, HunterLab;
6. Data: This mode can display all parameters that can be measured by the instrument except haze and hiding power;
7. Image: CIE LAB diagram, Yxy diagram, Luv diagram, reflection/transmittance diagram, K/S curve diagram, absorbance curve diagram;
9. Metamerism: measure metamerism parameters;
10. Liquid color: measure saybolt, ASTM color, platinum-cobalt color, Gander color;

The screenshot shows a software interface for haze measurement. The main window has a top bar with a home icon, 'Setting', and 'ASTM C/2°'. Below this are tabs for 'Instrument', 'Parameter', 'Tolerance', 'Display', 'Other', 'Guide', and 'Apply'. A left sidebar lists measurement types: 'H/T/C', 'Haze/T' (highlighted), 'Haze', 'Transmittance', 'Clarity', 'Color Diff.', 'Data', 'Figure', 'Metamerism', and 'Liquid'. The main area contains an inset window titled '标样测量 - 标样0098' and '试样测量'. This inset window features a graph of '透过率曲线(%)' with a y-axis from 0 to 200 and an x-axis for '波长 (nm)' from 400 to 700. A horizontal line is drawn at 100% transmittance. To the right of the graph is a table:

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On the right side of the inset window, there are icons for '设置' (Settings), '保存' (Save), and '测量' (Measure).

F-6

Other settings

1. The average setting window can set single measurement or average measurement;
2. The save setting window can be set to save manually or automatically;
3. The naming rule window can set the naming rule when the standard sample is saved.

The screenshot displays the 'Setting' application interface. At the top, there is a home icon, the title 'Setting', and 'ASTM C/2°'. Below the title bar are tabs for 'Instrument', 'Parameter', 'Tolerance', 'Display', 'Other', 'Guide', and 'Apply'. The 'Other' tab is active. It contains three main sections:

- Average:** Two radio buttons: 'Single Test' (selected) and 'Average Test'.
- Save Setting:** Two radio buttons: 'Manual Save' (selected) and 'Auto Save'.
- Naming Rules:** A table with two rows and two columns. The first row is for 'Target' and the second for 'Sample'. Each row has a text input field and a set of checkboxes for 'Number' (checked) and 'Date' (unchecked).

F-6

[Data browsing]

G

1. The left side of the page shows the standard sample data list, and the right side shows the sample data list under the standard sample;
2. There is a search box on the left side of the page, click to search the target data;
3. After clicking one of the standard samples, you can see the detailed information of the sample data under the standard sample data on the right side of the interface;
4. Long press the standard sample or sample to choose to recall, modify, delete the current selection, delete all;
5. Click the standard sample to enter the sample detailed information interface, you can search for the sample under the current standard sample, and you can export the current display data;
6. Click Parameter Edit to pop up the parameter editing window, where you can select the parameters displayed on the data interface.

Standard search box

Standard data list

Target Search

Target0001

Edit	Name	L*	a*	b*	dE*ab	P/F dB
Target	Target0001	99.99	-0.09	-0.06	--	--

Test Search Name ▾ Search Export Export All

G-1

- Target Search
- Target0002
- Target0001

Edit	Name	L*	a*	b*	dE*ab	P/F dB
Target	Target0002	99.99	-0.06	-0.08	--	--
0	Sample0001	100.01	-0.05	-0.13	0.05	Pas

Sample data list

Can export the current display data and export all data

Test Search Name Search Export Export All

G-2

Home Data View ASTM C/2°

Target Search

Target0002

Target0001

Edit	Name	L*	a*	b*	dE*ab	P/F dE
Target	Target0002	99.99	-0.06	-0.08	--	--
0	Sample0001	100.01	-0.05	-0.13	0.05	Pas

Set Into Target

Revise

Delete Selected

Delete All

Test Search Name Search Export Export All

G-3

ASTM C/2°

Data View

Preferences

Color space value	L*	Selected parameters	ADD DELETE REMOVE ALL	L* a* b* dE*ab P/F dE*ab	TOP UP DOWN BOTTOM
Color space	a*				
Color difference	b*				
Whiteness	c*				
Yellowness	h				
Blackness	X				
Transmittance	Y				
Color fastness	Z				
Force	x				
Color density					

FINISH

G-4

[Personal center]

H-1

The personal center interface can modify the account password and cancel the current account.

The screenshot displays a web interface titled "Individual Center" with a home icon and "ASTM C/2°" in the top right. The interface is divided into two main sections:

- Account: admin**: This section contains four input fields for "Company Name", "Address", "Contact Name", and "Email".
- Change Password**: This section contains three input fields for "Old Password", "New Password", and "New Password Confirm", followed by an "Enter" button.


At the bottom right of the interface, there is a "Logout" button.

H-1

[About]

I-1

About the interface, you can view the instrument information, such as software version, instrument version, instrument serial number, instrument model, etc.


I-1

[Log]

J-1

On the log interface, you can see the instrument login information, calibration information, instrument error information, etc.

The screenshot shows a web interface titled 'Daily Record' with a home icon on the left and 'ASTM C/2°' on the right. The main content area contains a list of log entries:


2016-01-02 20:43:05.3	100% successful calibration of instrument
2016-01-02 20:42:14.7	0% calibration succeed
2016-01-02 20:39:30.2	Account Login admin
2016-01-02 20:38:27.4	Account Login admin
2016-01-02 20:36:12.8	100% successful calibration of instrument
2016-01-02 20:35:56.4	0% calibration succeed
2016-01-02 07:34:28.4	100% successful calibration of instrument
2016-01-02 07:34:16.1	0% calibration succeed
2016-01-02 07:30:36.6	100% successful calibration of instrument
2016-01-02 07:30:23.9	0% calibration succeed

J-1

[Update]

K-1

In the case of networking, you can click Update to detect whether there is new software to obtain the latest software.


K-1

Introduction to the measurement interface

[Haze/transmittance/Clarity]

L-1

In this interface, the haze, total transmittance, sharpness and transmittance curve (%) can be displayed at the same time, and the difference between haze, total transmittance and sharpness can be calculated by comparing the standard samples, and passed simultaneously. The set tolerance automatically judges whether the sample is qualified.


L-1

[Haze/transmittance]

M-1

In this interface, the haze, total transmittance and transmittance curve (%) can be displayed at the same time, the difference between haze and total transmittance is calculated by comparing the standard sample, and the tolerance is automatically set at the same time Determine whether the sample is qualified.


M-1

[Haze]

M-1


In this interface, only the haze and transmittance curves (%) are displayed, and the difference in haze is calculated by comparing the standard sample, and at the same time, the sample is automatically judged by the set tolerance.

[Transmittance]

N-1

On this interface, only the total transmittance and transmittance curve (%) are displayed. The difference in total transmittance is calculated by comparing the standard sample, and the sample is automatically judged by the set tolerance.


N-1

[Clarity]

0-1

In this interface, only the clarity is displayed, the difference in sharpness is calculated by comparing the standard sample, and the sample is automatically judged by the set tolerance.

	Clarity
Target	100.00
Sample	100.00
Difference	0.00
Dudge	Pass

0-1

[Color difference]

P-1

CIELABCH

In this interface, the L^* , a^* , b^* , c^* , h values of the sample color can be measured, and the dL^* , da^* , db^* , dc^* , dH^* , dE^*_{ab} are calculated by comparing with the standard sample. At the same time, automatically determine whether the sample is qualified by the set tolerance.

Test Target --- Target0003

Test Sample --- Sample0001

Target

L* = 18.87
a* = 18.46
b* = 18.03
c* = 25.81
h = 44.34

Sample

L* = 19.03
a* = 18.56
b* = 17.25
c* = 25.34
h = 42.91

dL* = 0.15 Pass
da* = 0.10 Pass
db* = -0.78 Pass
dc* = -0.47 Pass
dH* = -0.64 Pass

dE*ab
0.81 Pass

Setting

Save


Measure

P-1

P-2

CIEDE2000

In this interface, the L^* , a^* , b^* , c^* , h values of the sample color can be measured, and the dL' , dC' , dH' , and dE^*2000 can be displayed by comparing with the standard sample calculation, and at the same time through the setting The tolerance of automatically determines whether the sample is qualified.

P-2

P-3

CIE94

In this interface, the L^* , a^* , b^* , c^* , h values of the sample color can be measured, and the dL^* , da^* , db^* , dc^* , dH^* , and dE^* are calculated by comparing with the standard sample. 94. At the same time, it automatically judges whether the sample is qualified by the set tolerance.

The screenshot shows a software interface for color measurement. At the top, it indicates 'ASTM C/2°'. Below that, there are two tabs: 'Test Target --- Target0003' and 'Test Sample --- Sample0001'. The interface is divided into three main sections:

- Target Data:**
 - $L^* = 18.10$
 - $a^* = 17.98$
 - $b^* = 16.59$
 - $c^* = 24.47$
 - $h = 42.70$
- Sample Data:**
 - $L^* = 17.79$
 - $a^* = 17.66$
 - $b^* = 15.99$
 - $c^* = 23.82$
 - $h = 42.16$
- Difference Values:**
 - $dL^* = -0.31$ Pass
 - $da^* = -0.32$ Pass
 - $db^* = -0.60$ Pass
 - $dc^* = -0.64$ Pass
 - $dH^* = -0.23$ Pass
- Final Result:**
 - $dE^*_{94} = 0.46$ Pass

On the right side of the interface, there are three buttons: 'Setting' (gear icon), 'Save' (floppy disk icon), and 'Measure' (crosshair icon).

P-4

CMC

In this interface, the L^* , a^* , b^* , c^* , h values of the sample color can be measured, and the dL^* , da^* , db^* , dc^* , dH^* , and $dE_{cmc}(l:c)$, at the same time automatically determine whether the sample is qualified by the set tolerance.

The screenshot displays a software interface for color measurement. At the top, it shows 'ASTM C/2°'. Below this, there are two tabs: 'Test Target --- Target0003' and 'Test Sample --- Sample0001'. The interface is divided into three main sections:

- Target Data:**

L^*	=	18.33
a^*	=	18.35
b^*	=	16.26
c^*	=	24.52
h	=	41.55
- Sample Data:**

L^*	=	18.32
a^*	=	18.51
b^*	=	16.65
c^*	=	24.90
h	=	41.97
- Difference Values:**

dL^*	=	-0.01	Pass
da^*	=	0.16	Pass
db^*	=	0.39	Pass
dc^*	=	0.38	Pass
dH^*	=	0.18	Pass

At the bottom, a box displays the final result: $dE_{cmc}(2.0:1.0)$ 0.30 Pass. On the right side, there is a vertical toolbar with icons for 'Setting', 'Save', and 'Measure'.

P-4

Hunter Lab

In this interface, the values of Hunter L, Hunter a, and Hunter b of the sample color can be measured, and dHunter L, dHunter a, dHunter b, and dEab can be calculated by comparing with the standard sample. At the same time, whether the data is automatically determined by the set tolerance qualified.

The screenshot displays the Hunter Lab software interface. At the top, it shows 'ASTM C/2°' and a home icon. Below this, there are two tabs: 'Test Target --- Target0003' and 'Test Sample --- Sample0001'. The main area is divided into three sections:

- Target (Hunter):**
 - L = 15.99
 - a = 11.51
 - b = 7.24
- Sample (Hunter):**
 - L = 15.92
 - a = 11.51
 - b = 7.23
- Comparison Results:**
 - dL = -0.07 **Pass**
 - da = -0.00 **Pass**
 - db = -0.00 **Pass**
 - dEab **0.07 Pass**

On the right side, there is a vertical toolbar with three icons: a gear for 'Setting', a floppy disk for 'Save', and a crosshair for 'Measure'.

[Data]

Q-1

1. In the data interface, you can click parameter edit to select the parameter you want to see;
2. Check the parameter difference of the sample by measuring the standard sample and then measuring the sample;
3. Click on the data to select it, and long press on the data to delete and rename the data.


Edit	Name	L*	a*	b*	dE*ab	P/F dE*ab
Target	Target0003	19.20	18.68	17.05	--	--
1	Sample0001	18.88	18.68	17.27	0.39	Pass

[Figure]

R-1

CIELAB


In this interface, you can measure the L^* , a^* , and b^* values of the sample, and at the same time use the a^* and b^* values of the sample to trace points on the CIELAB graph and display the L^* , a^* , and b^* values of the data.

R-2

Yxy


In this interface, you can measure the Y, x, and y values of the sample, and use the x and y values of the sample to trace points on the Yxy graph and display the Y, x, and y values of the data.

R-3

Luv


In this interface, you can measure the L^* , u^* , and v^* values of the sample, and use the u' and v' values of the sample to trace points on the Luv chart and display the L^* , u^* , and v^* values of the data.

R-4

k/s curve

In this interface, the k/s value of the sample can be measured, and the K/S curve under 360-780nm is displayed at the same time.


R-4

R-5

Transmittance curve


In this interface, the transmittance value of the sample can be measured, and the transmittance curve under 360-780nm can be displayed.

R-6

Absorbance curve


In this interface, the absorbance value of the sample can be measured, and at the same time, the absorbance curve under 360-780nm can be displayed.

[Metamerism]

S-1

The left side of the interface is the value calculated by the measurement sample using the first light source/angle, and the right side is the value calculated by the measurement sample using the second light source/angle. Metamerism values.

[Liquid color]

T-1

The cuvette optical path on the left of the interface is the recommended cuvette optical path size for measuring different parameters (for example, if you want to measure the saybolt parameter, the recommended cuvette optical path size at this time is 50mm), the right is the parameter value, and whether it is qualified Judgment.

Cuvette Light Path	Parameter	Target	Sample	Difference	Dudge
10mm	Pt-Co/Hazen/ APHA	500.00	500.00	0.00	Pass
10mm	Gardner Color	13.91	13.87	-0.04	Pass
50mm	Saybolt	-16	-16	0	Pass
33mm	ASTM Color	4.6	4.6	-0.0	Pass

T-1

Parameter introduction


Haze

Diffuse scattering reduces the image quality of the object. Small particles inside the material or sample surfaces cause scattering, and scattered light is scattered to different angles and the optical density at each angle is small. It causes a decrease in contrast and the sample forms a milky or cloud-like appearance, this phenomenon is known as haze.

Conditions for the assessment of transparency



Transparent product appearance has the characteristics of gloss, color and transparency. Transparency is particularly important, and its evaluation conditions are: transmittance, haze and so on. The ratio of light to incident light when the transparency rate is all projected. It will decrease as the surface of the material reflects and absorbs light. According to the ASTM D1003 the percentage of light that when passing through that deviates from the incident beam by greater than 2.5 degrees on average is defined as haze.

Measuring principle


When a beam of parallel light hits the sample and enters the integrating sphere, it is scattered on the sphere inner white coating, the total transmittance is measured when the optical trap outlet on the sphere is closed. When the optical trap is opened, haze will be measured.

Compensation method

The above picture is not using the compensation method to measure the transmission rate, in two times measurement, because the integral sphere's area is different, causes the integral sphere efficiency are different, thus has the certain influence to the measurement result accuracy.


The above picture is used to measure the transmission rate by the compensation method, a compensating port is added to the integral sphere, the first measurement sample is placed at the compensation port, the second measurement sample is placed at the measuring port, while the compensation port is in the open state, the integral sphere area of the two measuring process is identical, the integral sphere is consistent, so that the measurement result is more accurate.

Exception handling analysis


Error	Analysis	How to solve?
1.The instrument does not start.	Power connection may be abnormal	Check the Power interface for good contact and plug in the power supply.
2.No access to main interface after start.	The power-on calibration process may be abnormal	Re-calibrate as required
3.Error in measurement results.	Tolerance settings may be abnormal	Check tolerance settings and adjust
4.Test result is not correct	<ol style="list-style-type: none"> 1.The sample is close to the test port or not 2.Whether the sample surface is with scratches 	<ol style="list-style-type: none"> 1.Check the sample and test port fit to ensure close fitting 2.Check the sample surface condition to ensure that the sample is at good condition and has no effect on the measurement.

Annex


Standard accessories


Power Adapter


USB Cable

PC Software CD


Compensation Port Cover


0% Calibration Cover

Select Accessories

Film Sample Fixture

Liquid Sample Fixture

Company statement

- Our company commits to our customers 1 year warranty period for our Haze Meter series products from the date of the purchase, and our company shall be responsible to provide free maintenance for non-human caused malfunctions under normal usage. For malfunctions that are out of warranty period or caused by human factors, the company shall provide maintenance, and materials and repair shall be chargeable.
- The company is not liable for any loss or claim arising from the use of this product by the third party.
- The company is not liable for any damage caused by loss of data due to failure, maintenance or power off. To prevent the loss of critical data, be sure to back up all your data.
- The copyright of all products preset in this product belongs to the company and is protected by copyright law of the People's Republic of China.
- Our company's sale of this product does not represent the transfer or grant of any rights related to the copyright of the works, to the user.
- The product specifications and information mentioned in this specification are for reference only and will be updated at any time without prior notice.